




## *Consumables and Accessories*



# EDITORIAL

Since its creation in 1985, VLAD is an undeniable reference in autonomous energy : **Accumulators and Batteries for medical & industrial applications.**

We propose you a new range of medical Consumables and Accessories in addition to our cells & batteries offer. We invite you to have a look on this new catalogue that will guide you in making your choices.

Given the extent of our range, not all of our references are mentioned in this catalogue. **Do not hesitate to contact us; our team will be happy to answer and advise you according to your equipment and constraints.**

**Within 24 hours and in 95% of the cases, we can produce any specific battery, regardless of your field of activity.**








Our permanent stock of more 5,000 references over an inventory of 13,500 models allows us **to ship more than 90% of the orders we receive on the same day - for D+2 express deliveries anywhere around the world.**

We would like to thank you for your trust and interest in our catalogue.

The VLAD team

## OUR EXPERTISE



	<b>ECG CABLES</b>	<b>P.03</b>
	<b>SPO2 SENSORS</b>	<b>P.04</b>
	<b>OXYGEN CELLS</b>	<b>P.06</b>
	<b>MEDICAL LAMPS</b>	<b>P.09</b>
	<b>ECG PADS</b>	<b>P.12</b>
	<b>DEFIBRILLATOR PADS</b>	<b>P.13</b>
	<b>OTHER ACCESSORIES</b>	<b>P.15</b>

YOUR SUMMARY IN DETAIL



OEM reference	VLAD reference	Details	Brand	Devices	Connection
2320P-I	77700922	3-wired IEC Monobloc cable with clamps for monitor BP88s	COLIN MEDICAL		
2323P-I	77700970	3-wired IEC Monobloc cable with clamps for monitor SC5000	DRAGER	SC5000, SC6000, SC7000, SC8000, SC9000, Infinity Delta, Infinity Delta XL, Infinity Gamma	
2386P-I	77700937	3-wired IEC Monobloc cable with clamps for 2000	GE HEALTHCARE	Dash 1000, Dash 2000, Dash 2500, Dash 3000, Dash 4000, Dash 5000, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Responder 2000, Responder 3000, Corometrics120, GE® Vivid 7 Ultrasound	
2340RP-I	77700941	3-wired IEC Monobloc cable with clamps for Passport 2	MINDRAY		
2309P-I	77700967	3-wired IEC Monobloc cable with clamps for Life Scope-i	NIHON KOHDEN	BSM-2300 Life Scope, BSM-4100 Life Scope P, BSM-5100 Life Scope A, BSM-9510 Life Scope M, BSM-9800 Life Scope S	
2385P-I	77700931	3-wired IEC Monobloc cable with clamps for monitor Viridia	PHILIPS	Viridia, Ominicare CMA, Merlin, ML65A, ML66A, 43100A, 43110A, 43120A, 78352C, 78354C, 78834C, 862474 C3, Agilent, Agilent 50XM, CX50, Envisor, Heartstart, IE33, IntelliVue, M1001A, M1002A, M1165A, M1166A, M1167A, M1175A, M1177A, M1204A, M1205A	
2327P-I	77700940	3-wired IEC Monobloc cable with clamps for LP12	PHYSIOCONTROL	Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20	
23112P-I	77700953	3-wired IEC Monobloc cable with clamps for DG5000	SCHILLER	DG4000	
2396P-I	77700950	3-wired IEC Monobloc cable with clamps for 90492	SPACELABS MEDICAL	Ultraview 90492, 91220, 90496, 91496	
2353P-I	77700955	3-wired IEC Monobloc cable with clamps for Micropaq	WELCH ALLYN	PROPAQ LT	

SPO2 SENSORS


OXYGEN CELLS

MEDICAL LAMPS

ECG PADS

DEFIBRILLATOR PADS

OTHER ACCESSORIES



## Available finishes

- 3 / 4 / 5 / 6 / 10 strands
- Pliers / Pressions / Banana plug
- AHA cable on quotation



TO ORDER SEPARATE TRUNK OR THE EXTENSION CORD

CALL US FOR MORE INFORMATION

**REUSABLE - ONE PIECE - ADULT RIGID - DIGITAL**

Available in disposable x24


OEM reference	VLAD reference	Brand	Compatible devices	Technology	Sensor length	Connection
<b>U410-15R - M - ZOLL</b>	77700361	ZOLL	M series AED	MASSIMO	3m	
<b>U410-03 - M - SPACELAB</b>	77700412	SPACELABS MEDICAL	90419/ 01, 90651A/ 07, 90501, Ultraview 90496	NOVAMETRIX	3m	
<b>U410-15R - M - SCHILLER</b>	77700370	SCHILLER	Propaq CS, Encore, Micropaq (New) (M), ARGUS OXM PLUS ARGUS® VS-100	MASSIMO	3m	
<b>U410-15R - M - PHYSIOCONTROL</b>	77700369	PHYSIOCONTROL	Lifepack 20 (M)	MASSIMO	3m	
<b>U410-20 - M - PHILIPS</b>	77700420	PHILIPS	78352A, 78352C, 78354A, 78354C, 78834C, M1020A, M1025B, M1025A M1350B, M1772A/B, M1732A/B, M2475B	NELLCOR / NON OXIMAX	3m	
<b>U403-01 - M - NE TECH</b>	77700022	NELLCOR	NPB-40, NPB-20, NPB-75, NPB-180, NPB-185, NPB-190, NPB-195, NPB-290, NPB-295, NPB-3900, NPB-4000, NPB-3910, NPB-3920, NPB-3930, NPB-3940, N-10, N-20, N-30, N-100, N-180, N-190, N-200, N-3000, N-390, N-395, N-6000	NELLCOR / NON OXIMAX	1.1m	
<b>U410-48P - M - MINDRAY</b>	77700819	MINDRAY	Beneview T5/T8 Benevision, IPM/IMEC, VS-600/900	NELLCOR / OXIMAX	3m	
<b>U410-91 - M - GOLDWAY</b>	77700394	GOLDWAY MEDICAL	UT6000A G3, G60	PHILIPS	3m	
<b>U410-21MR - M - GE HEALTHCARE</b>	77700382	GE HEALTHCARE	Dash 3000/ 4000/ 5000, Solar 8000/8000M, Solar 8000i/9500 (Tram 451M, 851M), Cardiolab/Mac-Lab Corometrics 259, B20i, B40i, B60i Carescape B650, Carescape B850, Carescape V100	MASSIMO	3m	
<b>U403-01 - M - EDAN</b>	77700020	EDAN	H100B	NELLCOR / NON OXIMAX	1.1m	

## REUSABLE - ONE PIECE - ADULT RIGID - DIGITAL



Available in disposable x24

OEM reference	VLAD reference	Brand	Compatible devices	Technology	Sensor length	Connection
U403-01P - NM - DRAGER	77700047	DRAGER	Sirecust 1281, Monitor SC 6002XL, SC7000, SC9000XL	NELLCOR / OXIMAX	1.1m	
U403-06 - M - CORPULS	77700004	CORPULS	08/16	BCI	1.1m	
U410-46D - M - CONTEC	77700380	CONTEC	CMS60C/CW/D/DW, CMS60, PM50, PM60D CSM5000, CMS5000C, CMS5100, CMS6000, CMS7000/8000/9000, CMS5000, VET/6000-VET, CMS8000-VET/60D-VET		3m	
U403-01 - NM - COLIN	77700013	COLIN MEDICAL	BP306, BP88, BP508, BP808, M570		1.1m	
U410-15R - M - CAS MEDICAL	77700362	CAS MEDICAL	740-750	MASSIMO	3m	
U403-01 - NM - BIONET	77700011	BIONET	BM1, BM3, BM5 (New), BM3Plus/ BM3Wide		1.1m	
U410-80 - M - ARTEMA	77700353	ARTEMA	Diascope 8000, Athena, Diascope 1, Diascope 2, Diascope 3, Diascope G2	NON OXIMAX	3m	



## Other versions

- Adult / Pediatric / Neonate
- Finger / Wrap / Multi-site
- Flexible or Rigid



Available extensions



**OXYGEN CELLS**



OEM reference	VLAD reference	Product	Equivalence	Compatible devices
GO-06	77700816	Oxygen sensor GO-06	<b>ANALYTICAL</b> PSR-11-33-1, PSR-11-33-2 <b>ENVITEC</b> OOM107-2 <b>INTEGRAL</b> PROCESS 65000 <b>MAXTEC</b> MAX-18 (R116P40) <b>MERCURY</b> 10-103-01	
GO-07	77700811	Oxygen sensor GO-07	<b>ANALYTICAL</b> PSR-11-77 <b>DATEX</b> 6051-0000-219 <b>HAMILTON</b> HM-03 <b>INTEGRAL PROCESS</b> 65013 <b>MERCURY</b> 10-103-02 <b>NEWPORT</b> OM 100, VM 200	
GO-07S	77700817	Oxygen sensor GO-07s	<b>ANALYTICAL</b> PSR-11-60-03 <b>MAXTEC</b> MAX-3 (R109P11)	
GO-09	77700812	Oxygen sensor GO-09	<b>ANALYTICAL</b> PSR-11-37-D4 <b>DEVILBISS</b> O2ANA-101 <b>ENVITEC</b> OOM202-2 <b>MAXTEC</b> MAX-14 (R116P80) <b>MERCURY</b> 10-103-13	
GO-10	77700813	Oxygen sensor GO-10	<b>ANALYTICAL</b> PSR-11-915-4 <b>DATEX</b> 6050-0004-110 <b>ENVITEC</b> OOM110 <b>GE</b> M-10 (410014) <b>INTEGRAL</b> PROCESS 65008 <b>MAXTEC</b> MAX-10 (R112P10) <b>MERCURY</b> 10-103-10	<b>DATEX</b> 7900, Aespire, Aestiva 3000, Aestiva/5, Airsys, Avance, Engstrom Modulus SE7900, Excel 210-7900, Modulus SE1700, Smart Vent
GO-11	77700829	Oxygen sensor GO-11	<b>ANALYTICAL</b> PSR-11-917-2S <b>DATEX</b> 6051-0000-216	
GO-12	77700814	Oxygen sensor GO-12	<b>ANALYTICAL</b> PSR-11-917-MH1 <b>ATOM</b> C81828-R22M <b>DATEX</b> 0600-00-0002 <b>DRAGER</b> 8362030 E <b>ENVITEC</b> OOM102 <b>GE</b> R-22 M1223210-S <b>HAMILTON</b> HM-12 <b>HILLROM</b> R-22 <b>MAXTEC</b> MAX-12A, R109P48, MAX-23 (R116P06) <b>MERCURY</b> 10-103-15 <b>SPACELABS</b> 90514	<b>ATOM</b> V-2100G, V-2200 <b>DATEX</b> Centiva S/5 <b>DRAGER</b> Isolette 8000, Isolette C200, Isolette C2000 <b>GE</b> Carestation 30 <b>HAMILTON</b> Arabella, Multinex Analys, Raphael <b>HILLROM</b> C2000



OEM reference	VLAD reference	Product	Equivalence	Compatible devices
GO-13	77700831	Oxygen sensor GO-13	<b>ANALYTICAL</b> PSR-11-55 <b>DRAGER</b> 6736140 <b>INTEGRAL PROCESS</b> 65003 <b>MAXTEC</b> MAX-9 (R114R80)	
GO-14	77700830	Oxygen sensor GO-14	<b>ANALYTICAL</b> PSR-11-917-MHJ3 <b>MAXTEC</b> MAX-13C	
GO-15	77700827	Oxygen sensor GO-15	<b>ANALYTICAL</b> PSR-11-917-J <b>ATOM</b> OX-60 <b>DRAGER</b> MU16250 <b>ENVITEC</b> OOM102-1, OOM102-1-HEY <b>MAXTEC</b> MAX-7, MAX-13 (R115P10) <b>MERCURY</b> 10-103-18 <b>MINDRAY</b> 0600-00-0106 <b>SPACELABS</b> 015-0285-00 <b>AIR LIQUIDE</b> YR010100	<b>BIRD</b> 6400, 8400, VIP <b>DATEX</b> Engstrom Elvira, Engstrom Erica <b>MINDRAY</b> Anestar Plus <b>AIR LIQUIDE</b> Alys, Clarys, YR010100
GO-16	77700832	Oxygen sensor GO-16	<b>ANALYTICAL</b> PSR-11-917-M <b>CAREFUSION</b> 68289 <b>MAQUET</b> 6419332 <b>MAXTEC</b> MAX-16 (R114P70), MAX-16HT (R114P73) <b>MERCURY</b> 10-103-14 <b>PHILIPS</b> 1134668, 989805662341	<b>CAREFUSION</b> Avea, SIPAP <b>MAQUET</b> Servo 900, Servo 300 <b>RESPIRONICS</b> Trilogy EV300
GO-25	77700815	Oxygen sensor GO-25	<b>ANALYTICAL</b> PSR-11-917-J1 <b>MAXTEC</b> MAX-250ESF (R115P85) <b>MERCURY</b> 10-103-08	
GO-43	77700825	Oxygen sensor GO-43	<b>ANALYTICAL</b> PSR-11-915-G <b>MAXTEC</b> MAX-43 (R109P07) <b>MERCURY</b> 10-103-16 <b>DATEX</b> 6600-1278-600	<b>DATEX</b> Giraffe
GO-47	77700828	Oxygen sensor GO-47	<b>ANALYTICAL</b> RSR-11-58 <b>MAXTEC</b> MAX-47 <b>MERCURY</b> 10-103-20	<b>HUDSON</b> 5800, 5801
GO-48	77700834	Oxygen sensor GO-48	<b>ANALYTICAL</b> M-11-60, All-11-60 <b>MAXTEC</b> MAX-48	
MAX-250E	77700857	Oxygen sensor MAX-250E	<b>MAXTEC</b> MAX-250E, R125P03-002	<b>NEWPORT</b> E360, HT-70

ECG CABLES

SPO2 SENSORS

OXYGEN CELLS

MEDICAL LAMPS

ECG PADS

DEFIBRILLATOR PADS

OTHER ACCESSOIRES



OEM reference	VLAD reference	Product	Equivalence	Compatible devices
GO-01	77700807	Oxygen sensor GO-01	<b>ANALYTICAL</b> PSR-11-915 <b>DATEX</b> 0237-2034-700 <b>DRAGER</b> C-100 6803290 <b>ENVITEC</b> OOM101 <b>HAMILTON</b> HM-10, HM-13 <b>INTEGRAL PROCESS</b> 65001 <b>MAXTEC</b> MAX-1 (R113P90) <b>MERCURY</b> 10-103-17 <b>PHILIPS</b> 15201A	<b>DATEX</b> 125, 205, 51, 510, 5100, 5125, 5140, 5200, 5250, 5250 RGM, 5400, 7200, 7800, 7810, Excel 210-7800, Excel 7800, Modulus 4700, Modulus 7800, Modulus II 7810, Oxicap, S5172802
GO-02	77700808	Oxygen sensor GO-02	<b>ANALYTICAL</b> PSR-11-915-2, PSR-11-915-2i, G0-080 (Old CnS SKU) <b>DATEX</b> 6051-000-224 <b>DRAGER</b> 6803290, 6850645 <b>ENVITEC</b> OOM201, OOM104 <b>HAMILTON</b> HM-14, HM-13 <b>MAXTEC</b> MAX-11 (R113P10), MAX-11i (R113P06)	<b>DRAGER</b> 8000IC, I-8000 6000, Anemone, Babylog 8000, Evita 2, Evita 2 Dura, Evita 4, Evita Series, Fabius, Fabius GS, Fabius Series, Incubator 8000, Narkomed 1A, Narkomed 2, Narkomed 2B, Narkomed 3, Narkomed 4, Narkomed GS, Oxycom, Oxycom 100 D, Oxydig, PM8030, SA2
GO-03	77700833	Oxygen sensor GO-03	<b>ANALYTICAL</b> PSR-11-917-J6 <b>BIRD</b> 66030 <b>COVIDIEN</b> 4-020933-00 <b>DATEX</b> 6051-0000-217 <b>DRAGER</b> 6735142 <b>ENVITEC</b> OOM103-1, OOM202-1, OOM111 <b>INTEGRAL PROCESS</b> 65005 <b>INVIVO</b> 9445 <b>MAXTEC</b> MAX-13C (R115P01) <b>MERCURY</b> 10-103-05 <b>SMITHS</b> R-17	<b>BIRD</b> 6400, 8400, Sentry, VIP <b>COVIDIEN</b> 7820, PB 850 <b>DATEX</b> 7200, Engstrom Elvira, Engstrom Erica <b>RESPIRONICS</b> Esprit Lifecare
GO-04	77700809	Oxygen sensor GO-04	<b>ANALYTICAL</b> PSR-11-917-M2 <b>COVIDIEN</b> 4-072214-00 <b>DATEX</b> 6051-0000-222 <b>ENVITEC</b> OOM202 <b>HAMILTON</b> 396009, HM-12 <b>INTEGRAL PROCESS</b> 65016 <b>MAXTEC</b> MAX-12 (R109P45) <b>MERCURY</b> 10-103-06 <b>MAQUET</b> 9004979-E347E, 6419332-E380E <b>SPACELABS</b> 90513 AIR <b>LIQUIDE</b> YR049700	<b>BIRD</b> Avea <b>CAREFUSION</b> Viasys Avea 16448 <b>COVIDIEN</b> 740, 760, 840 <b>DRAGER</b> Infinity Delta <b>HAMILTON</b> 51254, Analyser, Multinex, Viasys <b>MAQUET</b> 9000, C/D, SV300, SV900, Servo 300, Servo 300A, Servo 700, Servo 710, Servo 900C, Servo 900D, Servo 900E, V710 <b>SPACELABS</b> Capnograph <b>AIR LIQUIDE</b> Cesar 2, Horus <b>VYAIR</b> Bellavista 1000
GO-05	77700810	Oxygen sensor GO-05	<b>ANALYTICAL</b> PSR-11-33 <b>HAMILTON</b> HM-01 <b>INTEGRAL PROCESS</b> 65007 <b>MAXTEC</b> MAX-19 (R116P60), R116P60 <b>MERCURY</b> 10-103-00 <b>AIR LIQUIDE</b> KV050501	<b>HAMILTON</b> Aladdin, Amadeus, Heyer Pocket 02, Veolar





AVAILABLE IN  
COMPATIBLE

MEDICAL LAMPS		SHADOWLESS					
OEM reference	VLAD reference	Device	Manufacturer	Voltage	Power	Base	Retail unit
CZ90522-0	77700783	C450	BERCHTOLD	22.8V	110W	G6.35	x 1
CZ90624-0	77700786	C570	BERCHTOLD	22.8V	150W	G6.35	x 1
CZ90722-0	77700788	C952	BERCHTOLD	22.8V	250W	G6.35	x 1
CZ90724-0	77700790	C950	BERCHTOLD	24V	250W	G6.35	x 1
CZ90822-0	77700792	D500	BERCHTOLD	22.8V	150W	G6.35	x 1
CZ90824-0	77700794	C571	BERCHTOLD	24V	150W	G6.35	x 1
CZ90922-0	77700796	D510	BERCHTOLD	22.8V	150W	G6.35	x 1
CZ95512-0	77700798	D300	BERCHTOLD	12V	50W	G6.35	x 1
56018550-0	77700767	Hamburg	HANAULUX	24V	40W	BA15D	x 1
56018566-0	77700769	London	HANAULUX	22.8V	50W	G6.35	x 1
56053010-0	77700771	Blue 30	HANAULUX	24V	50W	Special	x 1
56053026-0	77700773	Blue 80	HANAULUX	21.5V	130W	Special	x 1
56053198-0	77700775	Blue 90	HANAULUX	22.8V	90W	Special	x 1
56053296-0	77700777	Blue 100	HANAULUX	22.8V	113W	G6.35	x 1
56115282-0	77700779	2000	HANAULUX	22.8V	50W	G6.35	Pack x 8
567901947-0	77700781	3000	HANAULUX	22.8V	40W	G6.35	x 1
AX163508-0	77700765	ANGENIEUX	MAQUET	24V	150W	G6.35	x 1
AX186762-0	77700759	AXCEL	MAQUET	24V	100W	GX6.35	x 1
ECA002-0	77700764	PRISMATIC	MAQUET	24V	120W	X514	Pack x 10
LHX001C-0	77700761	PRISMALIX	MAQUET	23V	100W	X514SP	Pack x 10
8690301-0	77700800	ML301	MARTIN	12V	50W	G5.3	x 1
8690501-0	77700801	ML500	MARTIN	24V	120W	G6.35	Pack x 6
8690701-0	77700803	ML700	MARTIN	24V	150W	G6.35	Pack x 6
8691000-0	77700805	ML1000	MARTIN	24V	250W	G6.35	x 1

# LEADER IN EUROPE FOR CELLS & BATTERIES ORIGINALS / GENERICS



## BIOMEDICAL

- Defibrillators
- E.C.G.
- Patient lifts
- Surgical motors
- Respirators
- ...

## TECHNICAL

- Cells
- Accumulators
- SLA batteries
- Multimedia
- Fire-Rescue
- ...

Discover all  
our references on

[WWW.VLAD.FR](http://WWW.VLAD.FR)





MEDICAL LAMPS		OPHTHALMOSCOPY		
OEM reference	VLAD reference	Manufacturer	Voltage	Technology
X0288046	77700856	HEINE	3.5V	Xenon
X0288048	77700855	HEINE	3.5V	Xenon
X0288070	77700854	HEINE	3.5V	Xenon
X0288089	77700852	HEINE	3.5V	Xenon
X0288090	77700853	HEINE	3.5V	Xenon
WA03000	77700849	WELCH ALLYN	3.5V	Halogen
WA03100	77700909	WELCH ALLYN	3.5V	Halogen
WA03800	77700851	WELCH ALLYN	3.5V	Halogen
WA04900	77700850	WELCH ALLYN	3.5V	Halogen
WA08200	77700848	WELCH ALLYN	3.5V	Halogen

**+ 5 000**  
REFERENCES  
IN STOCK

Stock and prices at

**+33 (0)2 47 54 08 29**

or on


**www.vlad.fr**





ACCESSORIES & CONSUMABLES >

## ECG PADS



x 1

OEM reference	VLAD reference	Product	Brand	Dimensions	Gel	Support
<b>BLUESENSOR-R-SNAP</b>	77700990	Pack of 500 ECG pads BlueSensor R	AMBU	57 x 48	Liquid gel	Foam
<b>BLUESENSOR-VL-SNAP</b>	77700989	Pack of 500 ECG pads BlueSensor VL	AMBU	72 x 68	Liquid gel	Foam
<b>VLA-SC08</b>	77700995	Pack of 2 000 adults ECG pads	MULTIBrandS	50 x 50	Liquid gel	Foam
<b>VLA-SC40</b>	77700992	Pack of 2 000 adults ECG pads	MULTIBrandS	43 x 46	Liquid gel	Foam
<b>VLA-SF02</b>	77701010	Pack of 2 000 adults ECG pads	MULTIBrandS	41 x 46	Liquid gel	Foam
<b>VLA-SF13</b>	77700993	Pack of 2 000 adult/Pediatric ECG pads	MULTIBrandS	57 x 34	Liquid gel	Foam
<b>VLA-SM06</b>	77700994	Pack of 2 000 adults ECG pads	MULTIBrandS	53 x 53	Liquid gel	Microporous



# CREATE YOUR ACCOUNT

on our website [WWW.VLAD.FR](http://WWW.VLAD.FR)

- Up-to-date complete catalog
- Personalized price list
- Stock in real time
- Online Quote & Order



Available on  
**MOBILE**



Order by **WEB**  
**EDI** or **PUNCHOUT**



ACCESSORIES & CONSUMABLES >

**DEFIBRILLATOR PADS**



OEM reference	VLAD reference	Product	Manufacturer
9131-001	88881021	Adult pads for G3	CARDIAC SCIENCE
9730-002	88881022	Pediatric pads for G3	CARDIAC SCIENCE
XELAED001	88889242	Adult pads for G5	CARDIAC SCIENCE
XELAED003	88889243	Pediatric pads for G5	CARDIAC SCIENCE
CUA0512P	88881100	Pediatric pads for IPAD	CU MEDICAL
CUA0903PF	88881258	Adult pads for IPAD	CU MEDICAL
DDP-100	88881010	Adult pads for DCF200 / DCF210	DEFIBTECH
DDP-200P	88881019	Pediatric pads for DCF200 / DCF210	DEFIBTECH
989803149981	88881103	Adult pads pour FR3	PHILIPS
989803158211	88881040	Adult / Pediatric pads for FR2	PHILIPS
M5071A	88881038	Adult pads for HS1	PHILIPS



ECG CABLES

SPO2 SENSORS

OXYGEN CELLS

MEDICAL LAMPS

**ECG PADS**

**DEFIBRILLATOR PADS**

OTHER ACCESSORIES




Double



Unit or Pack x10

OEM reference	VLAD reference	Product	Manufacturer
M5072A	88881039	Neonate/Pediatric pads for HS1	PHILIPS
11101-000016	88889482	Pediatric pads for Lifepak 1000/500	PHYSIOCONTROL
11101-000021	88890006	Adult/Pediatric pads for CR2	PHYSIOCONTROL
11996-000017	88889481	Quick Combo adult pads for Lifepak 1000/500	PHYSIOCONTROL
96343	88881026	Adult/Pediatric pads for AED (before 2014)	PRIMEDIC
97085	88889811	Adult/Pediatric pads for AED (before 2014)	PRIMEDIC
SAV-C0016	88881037	Pediatric pads for SaverOne	SAVER ONE
SAV-C0846	88881036	Adult pads for SaverOne	SAVER ONE
0-21-0020	88881028	Adult pads for Fred Easy	SCHILLER
0-21-0037	88881029	Pediatric pads for Fred Easy	SCHILLER
0-21-0040	88890002	Adult pads for FRED PA-1	SCHILLER
0-21-0045	88890324	Pediatric pads for FRED PA-1	SCHILLER
8900-000260	88890004	Adult/Pediatric pads for AED3	ZOLL
8900-0801-01	88881032	Adult pads for AED+	ZOLL

**Wide choice available**

- Adult / Pediatric
- Original / Generic
- Pediatric / Child key

[➤ Consult us](#)



## MEDICAL BED ACCESSORIES



## REMOTE CONTROLS



OEM reference	Brand Device	VLAD reference	Functions	Device	Compatibility
8560026	ARJO	88889856	4 buttons = 2 functions	Alenti (Molex connector)	
8560023	ARJO	88889865	2 buttons = 1 function	Calypso	BOLERO / WALKER before 2009
TGO.310	ARJO	88889854	2 buttons = 1 function	Excelsior	
700-13700-box	ARJO	88890157	6 buttons / LCD Screen	MaxiMove (after 2007)	
8560020	ARJO	88889860	6 buttons = 3 functions	MaxiMove	
KTX00890.0	ARJO	88889853	2 buttons = 1 function	Minstrel	TRIXIE / SARITA
KTX11600.0	ARJO	88889857	4 buttons = 2 functions	Sara (SUB-D connector)	
8560004	ARJO	88889859	4 buttons = 2 functions	Alenti (hook connector)	



## MEDICAL BED ACCESSORIES



## CHARGERS



OEM reference	Brand Device	VLAD reference	Designation	Compatibility battery
NDA2200	ARJO	88881256	Wall charger for Alenti	NDA
NDA1000	ARJO	88889867	Wall charger for Ambulift	KKA
N/A	ARJO	88889870	Wall charger for Maxi-Twin	NEA
KPA8200	ARJO	88881144	Wall charger for Sara+	KPA
KTA0101	ARJO	88889810	Wall charger for Sarita	KTA
NDA8200	ARJO	88889419	Table charger for Alenti	NDA
N/A	ARJO	88889868	Table charger for Ambulift	KKA
700-24201	ARJO	88889871	Table charger for MaxiSky 440	MaxiSky 600
NEA1000	ARJO	88881351	Table charger for Maxi-Twin	NEA
N/A	ARJO	88889869	Table charger for Sara+	KPA
CHJ20000020A011	LINAK	88889928	Charger for Jumbo lifting system	BAJ1
1340100	MOLIFT	88889769	Charger for Powerpack	



OTHER CHARGERS AVAILABLE

CALL US FOR MORE INFORMATIONS

# OUR COMMITMENTS



## STATE-OF-THE-ART LOGISTICS

Thanks to our large stock and our state-of-the-art logistics, orders received before 4 pm can be processed for delivery the next day.



## THE BEST CHOICE

We are both a distributor of the major brands (Varta, Arts Energy, Saft, Panasonic, Yuasa, Philips, etc.) and a manufacturer of high-quality batteries. We have the product that best suits your needs.



## A STRONG COMMITMENT FOR THE ENVIRONMENT

Recycling used batteries, respect for the environment and controlling energy consumption have been at the heart of our environmental policy since 1985.



## ISO-CERTIFIED QUALITY

The quality of VLAD products & services has allowed us to obtain 3 ISO certifications for our quality management systems : ISO 9001, ISO 13485 (medical devices) and ISO 14001.

## CONTACT US

Tél. +33 (0)2 47 54 08 29

 [www.vlad.fr](http://www.vlad.fr)

